Sterile Implants & Instruments

Instructions for Use

Re+Line® Bunion Correction System

Please read this information carefully before using the Re+Line® Bunion Correction System.

DESCRIPTION

The Re+Line® Bunion Correction System consists of a sterile multi component bone fixation implant (Re+Line® Bunion Correction Plate and cortical screw). The Re+Line® Bunion Correction Plate is provided with a set of accessory instruments designed for preparation of the implant site, and insertion of devices into the bone. The cortical screw is packaged sterile but separate from the main system.

SYSTEM COMPONENTS

Cortical Headed Screw:	Surgical-grade implantable stainless steel
Re+Line® Bunion Correction Plate :	Surgical-grade implantable stainless steel
Drill Guide:	Surgical-grade stainless steel
Stabilizer (temporary) Pin: Drill Bit:	Surgical-grade implantable stainless steel
Handle/Driver:	Surgical-grade implantable stainless steel
	Polycarbonate and surgical-grade stainless steel

INDICATIONS

The Re+Line® Bunion Correction System is indicated for alignment, stabilization and fusion of fractures, osteotomies and arthrodesis of small bones such as the foot.

CONTRAINDICATIONS

- Patient conditions including insufficient quantity or quality of bone.
- Blood supply limitations and previous or active infections that may inhibit healing.
- Surgical procedures other than for the indications listed.
- Patients with conditions that limit their ability or willingness to follow postoperative care instructions.
- The device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature.
- Foreign-body sensitivity. Where material sensitivity is suspected, appropriate test should be made and sensitivity ruled out prior to implantation.

WARNINGS

This device should only be used by physicians with training in and a thorough understanding of orthopaedic surgery.

- The Re+Line® Bunion Correction System is for SINGLE USE ONLY. DO NOT RESTERILIZE OR REUSE
- The Re+Line® Bunion Correction System is not designed to withstand partial or full weight-bearing or excessive activity until healing has occurred. The post-op regimen prescribed by the physician should be strictly followed to avoid stresses applied to the implant.
- Do not use the implants or instruments appear damaged upon inspection
- Do not use in conjunction with implants of other metallic alloys
- Any decision to remove the device should take into consideration the potential risk to the patient.
- The Re+Line® Bunion Correction System has not been evaluated for safety and compatibility on the MR environment. The Re+Line® Bunion Correction System has not been tested for heating or migration in the MR environment.
- The Re+Line® Bunion Correction System must be used only with provided insertion instruments.

CAUTIONS

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. A detailed surgical technique in print and/or electronic formats can be obtained by contacting a Nextremity Solutions sales representative.
- Surgical instruments should only be used for their intended purpose.
- Leave adequate bone between the device legs and the osteotomy line

POTENTIAL ADVERSE EFFECTS

- Infection, deep and superficial
- Foreign body reactions
- Nonunion or delayed union which may lead to breakage of the implant
- Bending or fracture of implant
- Loosening or migration of the implant
- Nerve damage due to surgical trauma
- Inadequate healing
- Pain, soft tissue discomfort or abnormal sensation due to the presence of the device

PREPARATION PRIOR TO USE

- Inspect product for damage. Do not use if the product is damaged.
- Users of this system are encouraged to contact their Nextremity Solutions representative if, in their professional judgment, they require a more comprehensive explanation of the surgical technique to be used with this device.
- Insure that the patient is an appropriate candidate for the intended procedure based on the Indications, Contraindications, Warnings and Precautions listed above.
- Provide the patient with information concerning the need for limited weight bearing, and the possibility of device failure as the result of loosening, excessive activity or delayed or incomplete healing.

DIRECTIONS FOR USE

WARNING: Use only the insertion instruments provided.

- Sizing the Implant - Prior to surgery, assess pre-op x-ray to determine implant size. The size of the Re+Line® implant is determined by the width of the resected metatarsal head and the length of the flattened metatarsal shaft. The Re+Line® implant works most effectively when placed on a flat medial surface after correction and secondary trimming. The length of the medial flat surface will vary case-by-case, but generally the greater the correction of the capital fragment, the greater the need for a longer Re+Line® implant. As long as the metatarsal head medial lateral width is wider than 15 mm, a 20 or 25 mm implant should be used. The length of the resected metatarsal shaft determines between the 20 and 25 mm implant. The 20 and 25 mm Re+Line® Bunion Correction Plate are more commonly used.
- Prepare osteotomy site. Resect medial eminence mass as required. *Note: To reconfirm implant length, use a surgical metric from center of metatarsal head to desired proximal fixation point.*
- After assessing the size of the metatarsal head and deformity, open the appropriate size Re+Line® Kit. Center the Landmark® Positioning Guide on the medial side of the metatarsal head using the apex Landmark® (laser ring) as a guide.
- Pin the Landmark® Positioning Guide to the metatarsal, proximal to the osteotomy cuts. Drill the pin to stabilize the Landmark® Positioning Guide against the metatarsal. *Note: Fluoroscopy may be used to ensure placement. The Landmark® Positioning Guide functions as a template for the implant.*
- Using the 1.4 mm (shortest drill), drill the three remaining holes on the Landmark® Positioning Guide bushing. *Note: When drilling for the times, drill to the laser marking depth line. When drilling the apex hole, drill through the medial surface. Time hole depth is 15mm when the groove is at the top of the Landmark® Positioning Guide bushing.*
- (Optional) Mark the cutting landmark points on the metatarsal. Angle is 60 degrees between pen notches and apex.
- Remove the pin used to stabilize the Landmark® Positioning Guide. Cut the Chevron osteotomy and translate. Translate the head in a lateral direction to the desired position.
- Resect medial prominence off the metatarsal. The medial surface of the resected medial metatarsal head and metatarsal shaft should be level for proper seating of the Re+Line® Bunion Correction Plate. In order to prevent the head from moving, a distal to proximal K-Wire, inserted dorsally, may be used at this point to stabilize while cutting the redundant medial bone on the shaft.
- Insert the Re+Line® Bunion Correction Plate as far as possible into the pre-drilled holes on the medial surface of the metatarsal head. In hard bone, over-drilling the times may be necessary.
- While ensuring full reduction of the metatarsal head, maintain compression of the head at the osteotomy site on the metatarsal shaft. The Re+Line® Bunion Correction Plate should rest as flush against the metatarsal as possible. *Note: Fluoroscopy may be used prior to releasing plate to ensure proper placement.*
- Using the longer 2 mm drill, drill through the proximal end of the compression slot. Ensure the 2 mm pilot hole is made as proximal as possible to the center of the slot. This will ensure proper linear compression and thus reduction of the two bone fragments. With the drill tip crossing the lateral cortex, read the length marking that corresponds to the medial surface of the Re+Line® Bunion Correction Plate to select the appropriate screw length. Once accurate screw sizing is chosen, place the screw and tighten, reducing bone fragments.
- Using the 2 mm drill (longest drill), drill the second hole in the center; this is not for compression. After this drill hole is made and the proximal compression hole has been screwed and tightened, then insert second screw and tighten. Verify reduction, close in the usual manner.

Optional device removal: In the event the device needs to be removed after placement has been completed.

- Use the driver to carefully remove the screws from the bone.
- Upon removal of the screws gently remove the compression staple from the bone with the use of rongeurs.

CAUTION: IF THE DEVICE MUST BE REMOVED, DO NOT REUSE ANY PORTION OF THE DEVICE.

HOW SUPPLIED

The Re+Line® Bunion Correction System is available as a sterile surgical kit. This includes a compression staple and instruments, as well as a separate single sterile cortical headed screw packaged individually. The components of the Re+Line® Bunion Correction System are for single use only. **DO NOT RESTERILIZE. DO NOT REUSE. If product is re-used it may cause infection and may compromise the mechanical integrity of the product which may lead to failure.**

STORAGE AND HANDLING

Store in a cool, dry place and in a manner that protects the integrity of the packaging of all implants and instruments.

PACKAGING AND LABELING

- Nextremity Solutions devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.

Nextremity Solutions® and Re+Line® are trademarks of Nextremity Solutions, Inc. These products are covered by one or more U.S. and international patents pending. All rights reserved. Printed in the USA

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